



Product Catalog
2011

PREFACE

Dentapreg™ Products Catalog provides complete information about products offered by ADM, a.s.

Dentapreg™ products portfolio comprises following groups of products:

- Fiber reinforced composites (FRC) refills (page 3)
pre-impregnated light curing fiber reinforced composites.
- Accessories (page 11)
special hand instruments.

The products of each group are listed and described in separate chapter. The description covers basic information including end user price as well as picture of the product.

The **recommended end user prices** does not include VAT or/and any local taxes. Ordering codes are indicated by sign #.

Shipping terms are EXW.

Dentapreg™ products portfolio is a subject of the change, please check the current portfolio on www.dentapreg.com.

DENTAPREG™ FIBER REFILLS

Dentapreg™

Dentapreg™ Fiber Refills represent the core of Dentapreg™ products portfolio. The refills comprises of sticky pre-impregnated light curing glass fiber strips of different architecture (uni-directional, multi-directional, ultra-flat, enhanced pliability, enhanced strength). Dentapreg™ PINPost Refills comprises already cured flexible and elastic thin glass fiber reinforced pins.

List of the products:

Refills (Economy Packs)

- # D-001 Dentapreg™ SFU
- # D-002 Dentapreg™ SFM
- # D-003 Dentapreg™ PFU
- # D-004 Dentapreg™ PFM
- # D-005 Dentapreg™ UFM
- # D-006 Dentapreg™ PFUI
- # D-007 Dentapreg™ Starter Kit
- # D-008 Dentapreg™ Combi Praxis
- # D-009 Dentapreg™ Combi Labor
- # D-010 Dentapreg™ Combi

Minirefills

- # D-201 Dentapreg™ Minirefill SFU
- # D-202 Dentapreg™ Minirefill SFM
- # D-203 Dentapreg™ Minirefill PFU
- # D-204 Dentapreg™ Minirefill PFM
- # D-205 Dentapreg™ Minirefill UFM

PINPost

- # D-401 Dentapreg™ PINPost 0,35 mm Small Refill
- # D-402 Dentapreg™ PINPost 0,45 mm Small Refill
- # D-403 Dentapreg™ PINPost Large Refill
- # D-404 Dentapreg™ PINPost 0,45 mm Large Refill
- # D-405 Dentapreg™ PINPost 0,35 mm Large Refill

Dentapreg™ Fibers Refills

ITEM # D-001

Dentapreg™ SFU



Price 170.00 EUR

Pre-impregnated strips of 4.200 S glass fibers for perio, and trauma splinting, space maintainers and anatomical posts

- Blister of 8 strips 2 x 0,3 x 50 mm
- Light save box
- Instructions, virtual DVD

ITEM # D-002

Dentapreg™ SFM



Price 170.00 EUR

Pre-impregnated strips of 4.200 S glass fibers for perio, and trauma splinting, space maintainers and anatomical posts

- Blister of 8 strips 2 x 0,3 x 50 mm
- Light save box
- Instructions, virtual DVD

ITEM # D-003

Dentapreg™ PFU



Price 170.00 EUR

Pre-impregnated strips of 8.400 S glass fibers for frames of minimal invasive bridges, repair of removable dentures and anatomical posts.

- Blister of 8 strips 3 x 0,3 x 50 mm
- Light save box
- Instructions, virtual DVD

ITEM # D-004

Dentapreg™ PFM



Price 170.00 EUR

Pre-impregnated strips of 10.700 E glass fibers for frames of minimal invasive bridges.

- Blister of 8 strips 3 x 0,3 x 50 mm
- Light save box
- Instructions, virtual DVD

Dentapreg™ Fibers Refills

Dentapreg™

ITEM # D-005

Dentapreg™ UFM



Price 170.00 EUR

Pre-impregnated strips of 10.500 E glass fibers for strengthening crowns and large restorations and for repair and reinforcing of removable dentures.

- Blister of 8 strips 6 x 0,1 x 50 mm
- Light save box
- Instructions, virtual DVD

ITEM # D-006

Dentapreg™ PFUI



Price 255.00 EUR

Pre-impregnated flat, uni-directional strips of 8.400 S glass fibers for support of overdentures

- Blister of 4 strips 3 x 0,3 x 150 mm
- Light save box
- Instructions, virtual DVD

ITEM # D-007

Dentapreg™ Starter Kit



Price 190.00 EUR

Pre-impregnated strips for different applications.

- 4 strips of SFM (2 x 0,3 x 50 mm)
- 2 strips of PFU (3 x 0,3 x 50 mm)
- 2 strips of PFM (3 x 0,3 x 50 mm)
- 1 PINPost 0,45mm (Tube of 12 pins Ø 0,45 x 20 mm)
- Light save box
- Instructions, virtual DVD

ITEM # D-008

Dentapreg™ Combi Praxis



Price 170.00 EUR

Pre-impregnated strips for different applications.

- 4 strips of SFM (2 x 0,3 x 50 mm)
- 2 strips of PFM (3 x 0,3 x 50 mm)
- 2 strips of UFM (6 x 0,1 x 50 mm)
- Light save box
- Instructions, virtual DVD

Dentapreg™ Fibers Refills

ITEM # D-009

Dentapreg™ Combi Labor



Price 170.00 EUR

Pre-impregnated strips for different applications.

- 2 strips of SFU (2 x 0,3 x 50 mm)
- 4 strips of PFU (3 x 0,3 x 50 mm)
- 2 strips of UFM (6 x 0,1 x 50 mm)
- Light save box
- Instructions, virtual DVD

ITEM # D-010

Dentapreg™ COMBI



Price 170.00 EUR

Pre-impregnated strips for different applications.

- 2 strips of SFM (2 x 0,3 x 50 mm)
- 2 strips of PFU (3 x 0,3 x 50 mm)
- 2 strips of PFM (3 x 0,3 x 50 mm)
- 2 strips of UFM (6 x 0,1 x 50 mm)
- Light save box
- Instructions, virtual DVD

Dentapreg™ Fibers Minirefills

Dentapreg™

ITEM # D-201

Dentapreg™ Minirefill SFU




Price 90.00 EUR

Pre-impregnated strips of 4.200 S glass fibers

- Blister of 3 strips 2 x 0,3 x 50 mm
- Instructions

ITEM # D-202

Dentapreg™ Minirefill SFM




Price 90.00 EUR

Pre-impregnated strips of 7.300 E glass fibers

- Blister of 3 strips 2 x 0,3 x 50 mm
- Instructions

ITEM # D-203

Dentapreg™ Minirefill PFU




Price 90.00 EUR

Pre-impregnated strips of 8.400 S glass fibers

- Blister of 3 strips 3 x 0,3 x 50 mm
- Instructions

ITEM # D-204

Dentapreg™ Minirefill PFM




Price 90.00 EUR

Pre-impregnated strips of 10.700 E glass fibers

- Blister of 3 strips 3 x 0,3 x 50 mm
- Instructions

ITEM # D-205

Dentapreg™ Minirefill UFM



Price 90.00 EUR

Pre-impregnated strips of 10.500 E glass fibers

- Blister of 3 strips 6 x 0,1 x 50 mm
- Instructions

Dentapreg™ PINPost Refills

Dentapreg™

ITEM # D-401

Dentapreg™ PINPost 0,35 mm Small Refill



Price 20.00 EUR

Cured elastic and flexible pins of 960 S glass fibers for building anatomical stress free posts and retentions

- Tube of 12 pins Ø 0,35 x 20 mm
- Instructions

ITEM # D-402

Dentapreg™ PINPost 0,45 mm Small Refill



Price 20.00 EUR

Cured elastic and flexible pins of 1.440 S glass fibers for building anatomical stress free posts and retentions

- Tube of 12 pins Ø 0,45 x 20 mm
- Instructions

ITEM # D-403

Dentapreg™ PINPost Large Refill



Price 70.00 EUR

Cured elastic and flexible pins of 1.440 and 960 S glass fibers for building anatomical stress free posts and retentions

- 2x Tube of 12 pins Ø 0,35 x 20 mm
- 2x Tube of 12 pins Ø 0,45 x 20 mm
- Instructions

Dentapreg™ PINPost Refills

ITEM # D-404

Dentapreg™ PINPost 0,45 mm Large Refill



Cured elastic and flexible pins of 1.440 S glass fibers for building anatomical stress free posts and retentions

- 4x Tube of 12 pins Ø 0,45 x 20 mm
- Instructions

Price 70.00 EUR

ITEM # D-405

Dentapreg™ PINPost 0,35 mm Large Refill



Cured elastic and flexible pins of 960 S glass fibers for building anatomical stress free posts and retentions

- 4x Tube of 12 pins Ø 0,35 x 20 mm
- Instructions

Price 70.00 EUR

DENTAPREG™ ACCESSORIES

List of the products:

- # C-305 Dentapreg™ Shield
- # C-306 Dentapreg™ Fork

Dentapreg™ Accessories

ITEM # C-305

Dentapreg™ Shield



Special tool for protection of the rest of the Dentapreg™ strip during the light curing

Price 20.00 EUR

ITEM # C-306

Dentapreg™ Fork



Special tool for manipulation with Dentapreg™ strip

Price 20.00 EUR

Form 110421
All rights reserved © 2011
ADM, a.s., U Vodárny 2, 616 00 Brno, Czech Republic
www.dentapreg.com